



Service Bulletin for Spring Daylight Saving Time changes

In the upcoming transition to Daylight Savings Time, clocks will be moved forward one hour on all computers and servers. From a local time perspective, when this time change occurs the hour 2AM will become 3AM.

PeriCALM Tracings, being a UTC (universal time) based system, will continue to work normally during the time change, including the inbound and/or outbound EHR interfaces.

Outbound values sent over the EHR interface are sent with local time, so all values will be correctly shown with a gap of one hour between 2AM to 3AM on your EMR system.

PeriGen urges its users to verify any data entered and all alerts that may be displayed in your EMR system to avoid the possibility of features that may not function correctly for a brief period after the resetting of the clocks.

For any follow up questions regarding this process, please contact the PeriGen 24/7 support line at 866-321-6788, by mail support@perigen.com or through your account executive.